



# VKG Camera KX-Kymo

## Technology

- USB3 connection to PC / W11
- No extra processor unit needed
- Data flow controlled by our app

## Chips Info

- 1) Navigation: CMOS Sony Pregius, Global shutter
- 2) Kymography: CMOS Sony Pregius, 6500 fps
- Prism: Beam splitter 10/90



## Connectivity

- Intranet: local network
- Internet: Teams, Meet
- DICOM, PACS, GDT

## Common Features

- Save Images/Videos
- Simple Reporting in PDF
- Image printing

## Controls

- 3 programmable buttons
- touch screen
- footswitch

**Released for EU, CE according to the MDR  
registered in EUDAMED**